



## D6 Desktop Eye Tracking Solution


The D6 Optics is ASL's most robust desktop remote eye tracking solution. It is designed to quickly and accurately track the gaze position on all participants from infants to senior citizens. The stimulus can be presented on a single stationary surface, such as a computer or video monitor, or can be stationary real world objects. The system allows the participant approximately one square foot of head movement which eliminates the need for head restraint.

### There are many advantages to the D6 desktop remote solution:

- The only head free solution that uses bright pupil technology, giving superior capture and contrast.
- Quick and accurate head movement compensation.
- Constant visual feedback throughout entire tracking sessions.
- Automated procedures with manual overrides for challenging participants.
- Versatility in stimulus display devices— monitor, TV, projected image, even real world objects.
- Software Development Kit (SDK), which provides access to eye tracker controller port, serial out port and data files recorded by ASL interface program.
- Is displayed as a cursor or cross hairs superimposed on the image being viewed. A videotape or digital recording of this image can be created as a permanent record.

### Data Output

The D6 Optics are designed to measure a participant's eye line of gaze with respect to a single stationary surface in the environment. You can also configure the system to track participants viewing real objects. ASL solutions provide you with constant feedback on the eye tracking system and your participant's eye gaze. The gaze point is displayed as a cursor or cross hairs superimposed on the image being



Applied Science Laboratories  
175 Middlesex Turnpike, Bedford MA 01730 USA

Tel.: (781) 275-4000 Fax: (781) 275-3388 [www.asleyetracking.com](http://www.asleyetracking.com) Email: [asl@asleyetracking.com](mailto:asl@asleyetracking.com)



## D6 Desktop Eye Tracking Solution

A videotape or digital recording of this image can be created as a permanent record. Auto-calibration can be performed without operator intervention.

Recorded data includes time, x and y eye position coordinates, and pupil diameter. External data event markers can be recorded along with eye tracker data. Eye position coordinates correlate to specific areas on the surface being viewed. The D6 Optics includes a video head tracker, allowing head position to be recorded as well. The video head tracker does not require any equipment to be placed on the participant.

The D6 EYE-TRAC<sup>®</sup>6 operating software provides the system operator with the ability to enter calibration and participant data and specify the operating parameters of the D6 Optics. The provided software also converts data to text format for input to spreadsheets or other third party applications. Real-time digital data is available directly from the EYE-TRAC<sup>®</sup>6 control unit through a RS232 serial port.

### Optional Equipment

The D6 Optics is part of the ASL EYE-TRAC<sup>®</sup>6 Series. The ASL EYE-TRAC<sup>®</sup>6 Series offers the most eye tracking configuration options. This provides the greatest flexibility for your research requirements today and tomorrow.

- The EYE-TRAC<sup>®</sup>6 Control unit will work with all of the ASL EYE-TRAC<sup>®</sup>6 Optics options.
- New optics can be added at any time.
- System components and software can be shared among collaborators for greater allocation of time, funding, and resources.

D6 optics are portable and can be configured with a laptop as well as a PC.

### Training & Technical Support

ASL is committed to assisting researchers before, during and after the eye tracking data acquisition.

Unlimited technical support and free access to updates on the interface software are available at all times. Multiple licenses of the interface software are available at no additional charge.

We offer free scheduled training at our Bedford, MA (Boston) location for the life of the equipment. On site training is also available. Please check out our website for upcoming training sessions.

Technical Specifications	
Sampling Rate	60 Hz
Measurement Method	Pupil-corneal reflection
Accuracy	0.5° Visual angle
Resolution	0.25° Visual angle
Head Range	17" X 8" X 14"
Distance Range	20" to 40"
Dimensions	4.5" / 9.75" / 10.25"
H/W/D	11.43 / 24.765 / 26.035 cm
Visual Range	50° horizontal, 40° vertical
Real Time Data Outputs	X and Y gaze coordinates Horizontal and Vertical pupil diameter Two analog outputs One video output

Applied Science Laboratories  
175 Middlesex Turnpike, Bedford MA 01730 USA

Tel.: (781) 275-4000 Fax: (781) 275-3388 [www.asleyetracking.com](http://www.asleyetracking.com) Email: [asl@asleyetracking.com](mailto:asl@asleyetracking.com)